



**INTERNATIONAL STANDARD ISO 7197:2006**  
**TECHNICAL CORRIGENDUM 1**

Published 2007-07-01

INTERNATIONAL ORGANIZATION FOR STANDARDIZATION • МЕЖДУНАРОДНАЯ ОРГАНИЗАЦИЯ ПО СТАНДАРТИЗАЦИИ • ORGANISATION INTERNATIONALE DE NORMALISATION

## **Neurosurgical implants — Sterile, single-use hydrocephalus shunts and components**

### **TECHNICAL CORRIGENDUM 1**

*Implants neurochirurgicaux — Systèmes de dérivation et composants stériles, non réutilisables, pour hydrocéphalie*

#### **RECTIFICATIF TECHNIQUE 1**

Technical Corrigendum 1 to ISO 7197:2006 was prepared by Technical Committee ISO/TC 150, *Implants for surgery*, Subcommittee SC 3, *Neurosurgical implants*.

---

*Page 5*

**8.2 g)** the value of pressure in the text should read:

“... if applicable, include a graph of valve pressure over ambient pressure [between 0 kPa and 4,903 3 kPa (0 mm and 500 mm water column)] ...”